Densiron® 68

Densiron® Xtra

BG	1	IT	13
CZ	2	LT	14
DE	3	LV	15
DK	4	NL	16
EE	5	NO	17
EN	6	PL	18
ES	7	PT	19
FI	8	RO	20
FR	9	RU	21
GR	10	SE	22
HU	11	SI	23
IS	12	SK	24

Densiron® 68 / Densiron® Xtra

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Densiron® 68/Densiron® Xtra

Heavier-than-water endotamponades

Densiron® 68 / Densiron® Xtra are colourless homogeneous liquids composed of a mixture of ultrapure polydimethylsiloxane with the formula $(\text{CH}_3)_2\text{SiO}-[\text{Si}(\text{CH}_3)_2\text{O}]_n-\text{Si}(\text{CH}_3)_3$ and perfluorohexyloctane $(\text{C}_6\text{F}_{13}\text{C}_8\text{H}_{17})$. They are chemically and physiologically inert. Heavier-than-water endotamponades produced by Fluoron GmbH are characterised by their consistently high quality. Densiron® 68 / Densiron® Xtra are used as temporary retinal tamponades following successful retinal attachment and are removed dependent on the clinical findings. Densiron® 68 / Densiron® Xtra are used particularly in severe cases of retinal detachment for which conventional treatments would fail.

Indications

Densiron® 68 / Densiron® Xtra are used as temporary tamponades following surgical treatment of severe retinal detachments, especially in cases of:

- Inferior and posterior retinal holes
- Retinal detachments with giant tears
- Retinal detachments with proliferative vitreoretinopathy (PVR)
- Retinal detachments in cases of proliferative diabetic retinopathy (PDR)
- Traumatic retinal detachments



Composition

- 69.5 % polydimethylsiloxane
- 30.5 % perfluorohexyloctane

Physical properties

Device	Viscosity [mPas, 35 °C]	Density [kg/m³, 35 °C]	Refractive index [35 °C]
Densiron® 68	1200	1040	1.39
Densiron® Xtra	1000	1040	1.39

Spectral transmission

Operating instructions

Densiron® 68 / Densiron® Xtra are intended only for ocular application; they are applied at a dose determined by an ophthalmic surgeon for a defined indwelling time.

Preparation

Densiron® 68 / Densiron® Xtra in vials are drawn up undiluted into a sterile syringe with a suitable cannula attached.

Densiron® 68 / Densiron® Xtra in sterile glass syringes may be used immediately.

Attention: Densiron® 68 / Densiron® Xtra must be visually inspected for homogeneity and transparency prior to use. Products with phase separation or cloudiness must not be used!

Densiron® 68 / Densiron® Xtra are used as agents for temporary retinal tamponades following retinal reattachment. When removing Densiron® 68 / Densiron® Xtra from the eye it must be ensured that the heavier-than-water endotamponades are completely removed from the eye. Experience has shown that if the heavy silicone oil is removed in the post-operative phase, during which the visual acuity neither continues to increase nor drops again, the final result of the surgery can be improved. With the exception of eyes with a complex baseline situation, the indwelling time of Densiron® 68 / Densiron® Xtra must be kept as short as possible and restricted to a maximum of 3 months. Any silicone oil droplets of heavier-than-water endotamponades remaining in the eye or on an intraocular lens (IOL) can be dissolved with a semifluorinated alkane manufactured by Fluoron GmbH (e.g. F4H5® WashOut) and then removed.

Only medical personnel who have received training in the safe handling of the device are permitted to use the device.

Contraindications

Densiron® 68 / Densiron® Xtra are not permitted to be used in patients that are hypersensitive to silicone oil and perfluorohexyloctane or wear a silicone IOL. Because their density is greater than 1000 kg/m³, Densiron® 68 / Densiron® Xtra are not permitted to be used in cases of superior retinal detachment. In addition, laser coagulation and cryotherapy are not permitted on eyes filled with Densiron® 68 / Densiron® Xtra.

Side effects and complications

Densiron® 68 / Densiron® Xtra are generally well tolerated in the eye. The side effects and complications associated with a heavier-than-water endotamponade include optic nerve damage, corneal complications, secondary glaucoma, elevated intraocular pressure, hypotony, intraocular inflammations, emulsification, cataract, and decrease in visual acuity to vision loss.

Commercial form

Densiron® 68 / Densiron® Xtra are available in 10-mL units in vials and glass syringes.

Storage information

Densiron® 68 / Densiron® Xtra are not permitted to be stored below 2° C as this may cause phase separation, rendering the product unusable. Following storage for a week at room temperature, the phase separation will reverse and Densiron® 68 / Densiron® Xtra are then able to be used again.

Sterilisation

Sterility is ensured through dry heat and autoclaving. The device is autoclaved in the primary container (vial or glass syringe).

Warnings and precautions

- The indwelling time for Densiron® 68 / Densiron® Xtra should be kept as short as possible and limited to a maximum of 3 months.
- Before Densiron® 68 / Densiron® Xtra are permitted to be instilled in the eye, any heavy liquids, e.g. perfluorocarbons of any kind, must be completely removed from the vitreous chamber.
- Excessive instillation of Densiron® 68 / Densiron® Xtra in the posterior segment of the eye must be avoided.
- Glass syringes are potentially fragile due to their material structure. It is recommended to use appropriate application systems.
- Do not use this device after the expiry date. The expiry date is shown on the folding box.
- Vials and syringes are only intended for single use on one patient. Repeated use can lead to microbial contamination, which may cause a severe infection in the patient.

Please contact Fluoron GmbH if unexpected reactions or complications should occur during or following application that are suspected of being associated with Densiron® 68 / Densiron® Xtra.

Interactions

Interactions with Densiron® 68 / Densiron® Xtra may occur if perfluorocarbons are not completely removed (e.g. sticky oil). There are currently no known interactions with other substances.

Explanation of symbols

	Article number		Store dry
	Batch code		Sterilised with steam
	Manufacturer		Do not reuse
	Conformity mark and ID number of the Notified Body		Do not resterilise
	Follow instructions for use		Store between 10° C – 30° C
	Number of units		Crimped rim ampoule
	Use before		Syringe
	Do not use if packaging is damaged		Patient name or patient ID
			Healthcare institution/ provider

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FLUORON GmbH