Vioron®

BG	1	IT	13
CZ	2	LT	14
DE	3	LV	15
DK	4	NL	16
EE	5	NO	17
EN	6	PL	18
ES	7	PT	19
FI	8	RO	20
FR	9	RU	21
GR	10	SE	22
HU	11	SI	23
IS	12	SK	24

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Ladustada temperatuuril 0° C – 30° C

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Süstal

Vioron® is a biocompatible tissue dye for use as an aid in eye surgery. It is a sterile and pyrogen-free solution of trypan blue.

Indications

Vioron® was developed for ophthalmological surgical procedures involving the anterior segment of the eye (e.g. cataract operations and corneal grafts). The anterior lens capsule is made more readily visible through staining, thus simplifying capsulorhexis and minimising the risk of capsular tearing.

Vioron® also facilitates the dissection and transfer of the donor cornea in the case of lamellar corneal transplants and the removal of the abnormal Descemet membrane in the case of DMEK (Descemet Membrane Endothelial Keratoplasty) and DSA(A)EK [Descemet's Stripping (Automated) Endothelial Keratoplasty].

Composition


1 mL Vioron® contains:

- 0.6 mg trypan blue
- Disodium hydrogen phosphate ($\text{Na}_2\text{HPO}_4 \times 2 \text{H}_2\text{O}$)
- Monosodium phosphate ($\text{NaH}_2\text{PO}_4 \times 2 \text{H}_2\text{O}$)
- Sodium chloride (NaCl)
- Water for injections

Physical properties

Device	Density [kg/m³, 35°C]	Refractive index [35°C]
Vioron®	1000	1.33

Spectral transmission

Operating instructions

Vioron® is intended only for ocular application; it is applied at a dose determined by an ophthalmic surgeon for a defined indwelling time.

Preparation

Vioron® in vials is drawn up undiluted into a sterile syringe with a suitable cannula attached.

Suitable cannulas are attached directly to Vioron® in sterile glass syringes.

In the case of capsulorhexis, the intraocular fluid is first displaced from the anterior chamber with air or a viscoelastic substance in order to prevent Vioron® from becoming diluted. A cannula is then used to carefully apply Vioron® below the air bubble or viscoelastic substance to the anterior lens capsule within the ring of the iris.

Once adequate staining of the anterior lens capsule has been achieved, the anterior chamber is rinsed in order to remove any excess dye. The cataract operation can then proceed normally.

In the case of DMEK, the Descemet membrane containing the endothelium is able to be made more readily visible by staining it with Vioron® in order to remove the prepared donor cornea. The donor lamella obtained in this manner forms the Descemet roll, which is then able to be stained in vitro with Vioron®. Staining facilitates correct orientation of the roll, transfer to the recipient eye using an inserter, and subsequent unscrolling and application of the roll to the recipient stroma.

Exposure to the staining solution must be kept as short as possible.

Vioron® must be removed immediately following injection by completely flushing out the anterior chamber. The dissected Descemet roll must also be thoroughly rinsed following staining in vitro.

Only medical personnel who have received training in the safe handling of the device are permitted to use the device.

Contraindications

Although no allergic reactions to Vioron® have been reported, allergic reactions to the dye are unable to be ruled out in hypersensitive patients.

Use in pregnant women, nursing mothers, and children is not recommended.

Side effects and complications

Vioron® is generally well tolerated in the eye. If Vioron® is improperly injected into the anterior chamber, the dye may penetrate the vitreous chamber if the zonular fibres are relaxed and loose and temporarily stain the vitreous. In the case of capsulorhexis or corneal grafting (with application of vital dyes), corneal oedema, inflammatory reactions, increase or decrease in intraocular pressure, endothelial cell loss or clouding of the posterior lens capsule may occur postoperatively. No side effects/ complications associated with the dye have been reported to the manufacturer since the market launch of Vioron®.

Commercial form

Vioron® is available in 0.5-mL units in vials and glass syringes.

Storage information

Vioron® is stored at room temperature. Protect against light and freezing.

Sterilisation

Sterility is ensured through autoclaving. The device is autoclaved in the primary container (vial or glass syringe).

Warnings and precautions

- Do not use Vioron® after the expiry date. The expiry date is shown on the folding box.
- Use only undamaged packages containing clear, blue solution.
- Vials and glass syringes are only intended for single use on one patient. Repeated use can lead to microbial contamination, which may cause severe infection.
- Glass syringes are potentially fragile due to their material structure. It is recommended to use appropriate application systems.
- Prolonged exposure to Vioron® may damage the corneal endothelial cells. As a result, exposure must be kept short. The maximum application duration of Vioron® must be limited to 30 minutes. Tissue changes due to Vioron® are unable to be ruled out in the case of longer applications.
- Check to ensure that the syringe stopper moves smoothly prior to use. The stopper must be retracted prior to the injection.
- Use in pregnant women, nursing mothers, and children is not recommended.

Please contact Fluoron GmbH if unexpected reactions or complications should occur during or following application that are suspected of being associated with Vioron®.

Interactions

There are currently no known interactions with other substances.

Explanation of symbols:

	Article number		Do not use if packaging is damaged
	Batch code		Store dry
	Manufacturer		Keep away from sunlight
	Conformity mark and ID number of the Notified Body		Sterilised with steam
	Follow instructions for use		Do not reuse
	Number of units		Do not resterilise
	Use before		Store between 0° C – 30° C
			Crimped rim ampoule
			Syringe

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FLUORON GmbH