MD

Vioron®

BG	1	LT	14
CS	2	LV	15
DA	3	NL	16
DE	4	NO	17
EL	5	PL	18
EN	6	PT	19
ES	7	RO	20
ET	8	RU	21
FI	9	SK	22
FR	10	SL	23
HR	11	SV	24
HU	12	TR	25
IT	13		

FLU  RON®

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Vioron® is a biocompatible tissue dye for use in ophthalmic surgery. Vioron® is a sterile and endotoxin-free solution of trypan blue.

Composition


0.5 ml weakly buffered Vioron® solution [pH 7.0 - 7.8] contains:

- 0.30 mg trypan blue ($C_{36}H_{38}N_4Na_2O_{14}S_4$)
- 0.95 mg disodium hydrogen phosphate ($Na_2HPO_4 \times 2H_2O$)
- 0.15 mg sodium dihydrogen phosphate ($NaH_2PO_4 \times 2H_2O$)
- 4.10 mg sodium chloride (NaCl)
- water for injection

Physical properties

Device	Density [kg/m³, 25 °C ± 2 °C]	Osmolality [mOsm/kg]
Vioron®	990 – 1020	290 ± 15

Spectral transmission

Intended use

Vioron® is a sterile dye for selective staining of intraocular structures.

Indications

Vioron® is intended only for ocular application.

Vioron® was developed for ophthalmological surgical procedures involving the anterior segment of the eye (e. g., cataract operations and corneal grafts). The anterior lens capsule is made more readily visible through staining, thus simplifying capsulorhexis and minimising the risk of capsular tearing.

Furthermore, Vioron® is used to check the vitality of corneal grafts (corneal bank) for penetrating keratoplasty, as well as in lamellar keratoplasty (DMEK, DSAJEK) to improve the visibility of the donor preparation and to support the correct orientation of the endothelial Descemet membrane roll.

Expected clinical benefit

The performance characteristics of the Vioron® device are the following:

- improved visibility of the anterior lens capsule/endothelium descemet membrane (procedural success),
- successful/complete continuous curvilinear capsulorhexis/cornea transfer (anatomic success),
- improvement of visual acuity (functional success).

Use

Vioron® in sterile syringes is connected directly to a suitable cannula.

For capsulorhexis, the aqueous humor is first displaced from the anterior chamber with air or a viscoelastic to prevent dilution of Vioron®. Vioron® is then carefully applied with the cannula under the air bubble or the viscoelastic onto the anterior lens capsule within the iris ring. The indwelling time of Vioron® in the anterior chamber should be kept as short as possible and depends individually on the surgeon's decision. An indwelling time of seconds to about 1 minute is usually sufficient. The trypan blue solution is then rinsed out completely. Cataract surgery can then be performed as usual.

To check the vitality of corneal grafts for penetrating keratoplasty, Vioron® is used in the so-called "dye exclusion assay": living cells exclude the dye, dead cells are stained dark blue. In DMEK, Vioron® is used to improve the visualization of the thin transparent donor preparation.

Furthermore, staining with Vioron® supports the visualization and correct orientation of the endothelial Descemet membrane roll in the eye. Care should be taken to ensure the shortest possible exposure time to the staining solution.

Vioron® should be removed immediately after injection by complete rinsing of the anterior chamber. Likewise, after in vitro staining, thoroughly rinse the prepared Descemet roll.

Dosage

Vioron® is to be used undiluted. Vioron® in sterile syringes may be used directly.

Operating instructions

Vioron® is intended only for ocular application. Such application is permitted to be performed only by a trained ophthalmic surgeon.

Contraindications

Although no allergic reactions to Vioron® have been reported, allergic reactions to the dye cannot be ruled out in hypersensitive patients.

Side effects and complications

Vioron® is generally well tolerated in the eye. If Vioron® is injected improperly into the anterior chamber, the dye may penetrate into the vitreous cavity in the case of loose zonular fibers and temporarily stain the posterior lens capsule and vitreous. Side effects reported after the use of trypan blue include discoloration of intraocular lenses with high water content [hydrophilic acrylic lenses].

In capsulorhexis or keratoplasty, anterior eye segment complications, intraocular inflammations, decreased or increased intraocular pressure, endothelial cell loss, or posterior lens capsule opacification may occur postoperatively.

Warnings and precautions

- Do not use this device after the expiry date. The expiry date is shown on the devices and on the folding boxes.
- The device is only intended for single use on one patient. Repeated use may lead to microbial contamination, which may cause a severe infection in the patient.
- Use in pregnant women, nursing mothers and children is not recommended.
- Check the integrity of the sterile packaging. Do not use the device if the sterile packaging has been damaged.
- The indwelling time of the device should be kept as short as possible and limited to the recommended duration of use. A short length of indwelling time can reduce potential side effects and complications.
- A leak-free and safe connection of the syringe is ensured by the standardised Luer Lock [according to ISO 80369].
- Glass syringes are potentially fragile due to their material structure. It is recommended that appropriate application systems be used.
- Check that the syringe stopper moves smoothly prior to use. The stopper must be retracted prior to the injection
- If the exposure time is too long, Vioron® may damage corneal endothelial cells. Therefore, care should be taken to ensure the shortest possible exposure. In case of long-term application, tissue changes caused by Vioron® cannot be excluded.

Should unexpected side effects or complications occur during or following application and be suspected of being associated with Vioron, please contact Fluoron GmbH immediately via e-mail: complaints@fluoron.de or the competent authority of your country.

Interactions

There are currently no known interactions with other substances.

Sterilisation

Sterility is ensured by autoclaving. The device is autoclaved in the primary packaging.

Storage conditions

Vioron® is to be stored between 10 °C and 30 °C and is also to be protected from sunlight and humidity.

Commercial forms

Vioron® is available as a 0.5 ml single dose in suitable syringes.

Disposal

The device is to be disposed of in accordance with the regulations of your organisation. WHO recommendations: "Safe management of wastes from health-care activities, 2014, WHO"

Explanation of symbols

	Manufacturer		Store dry
	Use-by date		Store between 10 °C – 30 °C
	Batch number		Do not reuse
	Article number		Country symbol
	Distribution partner		Follow instructions for use https://ifu.fluoron.de/
	Sterilized using steam		Attention, important note
	Double sterile barrier system with protective packaging outside		Pyrogen-free
	Do not resterilize		Medical device
	Do not use if the sterile barrier system of the product or its packaging is damaged		Unique identifier of a medical device
	Fragile, handle with care		Quantity
	Protect from sunlight		Syringe
			Conformity mark

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FLUORON GmbH