MD

F-Octane

F-Decalin

BG	1	LT	14
CS	2	LV	15
DA	3	NL	16
DE	4	NO	17
EL	5	PL	18
EN	6	PT	19
ES	7	RO	20
ET	8	RU	21
FI	9	SK	22
FR	10	SL	23
HR	11	SV	24
HU	12	TR	25
IT	13		

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F-Octane and F-Decalin are colourless, homogenous liquids of high density composed of high purity perfluorocarbons. F-Octane and F-Decalin contain no biologically active components, are chemically and physiologically inert, non-toxic and laser stable.



Composition

F-Octane	F-Decalin
100 % perfluorooctane [C ₈ F ₁₈]	100 % perfluorodecalin [C ₁₀ F ₁₈]

Physical characteristics

Device	Density [kg/m ³ , 35 °C ± 2 °C]	Refractive index [35 °C ± 2 °C]
F-Octane	1730 ± 10	1.262 ± 0.001
F-Decalin	1910 ± 10	1.309 ± 0.001

Spectral transmission

Intended use

F-Octane/F-Decalin are sterile medical devices composed of ultrapure, liquid perfluorocarbon (perfluorooctane or perfluorodecalin). F-Octane/F-Decalin are used as intraoperative endotamponade in the field of ophthalmic surgery.

Indications

F-Octane/F-Decalin are intended only for ocular application.

F-Octane/F-Decalin are used for the treatment of severe cases of retinal detachment and for treatment of dislocated lenses, especially for:

- Retinal detachments (folded retina)
- Retinal detachments with giant tears
- Traumatic retinal detachments
- Lens dislocation

Expected clinical benefit

The performance characteristics of the F-Octane/F-Decalin devices are the following:

- easy and safe unfolding of detached retina, how well luxated lenses can be handled, tamponade can be removed well (procedural success),
- reattachment rate in cases of retinal detachment, cancellation of dislocation of the lens (anatomic success),
- improvement of visual acuity (functional success).

Use

After a partial or complete vitrectomy, F-Octane/F-Decalin are injected exactly opposite the optic disc. Due to their low viscosity, standard cannulas can be used. At the end of surgery, F-Octane/F-Decalin must be completely removed from the eye and, if necessary, replaced with a long-term tamponading agent.

Dosage

F-Octane/F-Decalin are to be used undiluted.

F-Octane/F-Decalin in vials are transferred undiluted into a sterile syringe with a suitable cannula attached.

F-Octane/F-Decalin in sterile syringes may be used directly.

Operating instructions

In the case of giant tears, prevent the fluid from penetrating the subretinal tissue.

F-Octane/F-Decalin are intended only for ocular application. Such application is permitted to be performed only by a trained ophthalmic surgeon who will determine the dosage and indwelling time.

Contraindications

F-Octane/F-Decalin must not be used in patients who are hypersensitive to perfluorocarbons. Due to their physical properties, F-Octane/F-Decalin are not permitted to be used as a long-term tamponade.

Side effects and complications

F-Octane/F-Decalin are generally well tolerated in the eye. The side effects and complications associated with F-Octane/F-Decalin include anterior eye segment complications, decreased or increased intraocular pressure, intraocular inflammations, emulsification, incomplete removal, formation of "sticky oil" in contact with silicone oil and subretinal penetration. Special caution should be taken when F-Octane/F-Decalin may enter the vascular system [e.g., during endoresection of intraocular tumors and ocular trauma]. Rare cases of gas embolism have been reported that may be associated with perfluorooctane/perfluorodecalin penetration.

Warnings and precautions

- Do not use this device after the expiry date. The expiry date is shown on the devices and on the folding boxes.
- The device is only intended for single use on one patient. Repeated use may lead to microbial contamination, which may cause a severe infection in the patient.
- Use in pregnant women, nursing mothers and children is not recommended.
- Check the integrity of the sterile packaging. Do not use the device if the sterile packaging has been damaged.
- The indwelling time of the device should be kept as short as possible. A short length of indwelling time can reduce potential side effects and complications.
- A leak-free and safe connection of the syringe is ensured by the standardised Luer Lock (according to ISO 80369).
- Glass syringes are potentially fragile due to their material structure. It is recommended that appropriate application systems be used.
- Overfilling of the posterior segment of the eye with F-Octane/F-Decalin should be avoided.
- Prior to unfolding the retina using F-Octane/F-Decalin, all traction must be removed in order to prevent the additional weight of the perfluorocarbons from inducing iatrogenic retinal detachments or retinal breaks.
- F-Octane/F-Decalin must not be injected directly into the vena vorticiosa. In addition, avoid inadvertent leakage of intraocular perfluorooctane/perfluorodecalin into the vena vorticiosa.

Should unexpected side effects or complications occur during or following application and be suspected of being associated with F-Octane/F-Decalin, please contact Fluoron GmbH immediately via e-mail complaints@fluoron.de or the competent authority of your country.

Interactions

Interactions with silicone oils may occur if F-Octane/F-Decalin are not completely removed [e.g., sticky oil]. There are currently no known interactions with other substances.

Sterilisation

Sterility is ensured by autoclaving. The product is autoclaved in the primary packaging material syringe or vial.

Storage conditions

F-Octane/F-Decalin is to be stored between 10 °C and 30 °C and is also to be protected from humidity.

Commercial forms

F-Octane/F-Decalin are available as 5 ml or 7 ml single doses in suitable vials and syringes

Disposal

The device is to be disposed of in accordance with the regulations of your organisation. WHO recommendations: "Safe management of wastes from health-care activities, 2014, WHO".

Explanation of symbols

	Manufacturer		Store dry
	Use-by date		Store between 10 °C–30 °C
	Batch number		Do not reuse
	Article number		Country symbol
	Distribution partner		Follow instructions for use https://ifu.fluoron.de/
	Sterilized using steam		Attention, important note
	Double sterile barrier system with outer protective packaging		Pyrogen-free
	Do not resterilize		Medical device
	Do not use if the sterile barrier system of the product or its packaging is damaged		Unique identifier of a medical device
	Fragile, handle with care		Quantity
	Protect from sunlight		Syringe
			Crimp rim vial
			Conformity mark

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FLUORON GmbH